

PM 200



Model	PM 200
Total magnification	100x/500x
Eyepiece	WF 10x (ø18mm)
Adapter	1x C-mount adapter
Plan achromatic objectives long working distance	PL L 10x/0.25 WD=7.30mm
	PL L 50x/0.65 WD=0.50mm
Focusing System	Coaxial coarse/fine focus system
Mechanical tube length	160mm
Reflected illumination	Coaxial illumination
	0.5W LED lamp rechargeable with brightness control
Polarization	Push/pull polarizer
Analyzer	Push/pull analyzer
Contrast techniques	Brightfield: Yes
	Darkfield: Not available
	DIC: Not available
	POL Contrast: Yes
Packing	Portable aluminium case LxWxH 330x270x200mm
Height	200mm
Weight	2.5kg

Optional accessories

Model	Description
PM O-200-21	Long working distance objective 20x/0.40 WD=8.60mm
PM O-200-22	Long working distance objective 60x/0.75 WD=1.34mm
PM O-200-23	Long working distance objective 80x/0.80 WD=0.96mm
PM O-200-24	Long working distance objective 100x/0.85 WD=0.40mm
PM S-200-61	Magnetic stage

PM 200

Long working distance objective 20x/0.40 WD=8.60mm