

**ADDITIONAL
INSTRUMENTS**



Digital ultrasonic
cleaners

Model	Type	Capacity	Heating	Frequency	Ultrasonic power	Heating power	Power supply	Tank dimensions	Overall dimensions	Weight	Recommended for:
SONIC-D06	Digital ultrasonic cleaner	0,6 l	No	40 kHz	50 W	-	230V - 50Hz	150x85x65 mm	175x110x155 mm	1,5kg	For cleaning samples surfaces after the cutting and mounting processes, before the grinding/polishing processes
SONIC-D06H	Digital ultrasonic cleaner	0,6 l	Yes	40 kHz	50 W	30 W	230V - 50Hz	150x85x65 mm	175x110x155 mm	1,5kg	
SONIC-D16H	Digital ultrasonic cleaner	1,3 l	Yes	40 kHz	70 W	30 W	230V - 50Hz	150x135x100 mm	180x165x205 mm	2,5kg	
SONIC-D32H	Digital ultrasonic cleaner	3,2 l	Yes	40 kHz	120 W	100 W	230V - 50Hz	240x135x100 mm	265x165x235 mm	4kg	
SONIC-D55H	Digital ultrasonic cleaner	4,5 l	Yes	40 kHz	180 W	200 W	230V - 50Hz	300x150x100 mm	330x180x260 mm	5,5kg	
SONIC-D70H	Digital ultrasonic cleaner	7 l	Yes	40 kHz	240 W	500 W	230V - 50Hz	300x240x100 mm	330x270x260 mm	8,5kg	
SONIC-D150H	Digital ultrasonic cleaner	15 l	Yes	40 kHz	360 W	500 W	230V - 50Hz	330x300x150 mm	360x330x330 mm	11kg	