

## ZETAPRESS M2D

## Features

- > Automatic mounting press with large size LCD touch screen, microprocessor controlled for setting all process parameters such as, automatic pre-heating, heating, working pressure, cooling time
- > One cylinder with intermediate piston which makes possible to carry out 2 mounting operations at same time
- > Different mold sizes available ( $\varnothing$  25mm,  $\varnothing$  32mm,  $\varnothing$  40mm,  $\varnothing$  50mm), two sizes can be selected
- > Automatic water cooling system

TWO DIFFERENT  
CHAMBERS DIAMETER



Model	ZETAPRESS M2D
Mold diameter (2 sizes can be selected as standard configuration)	$\varnothing$ 25, $\varnothing$ 30, $\varnothing$ 40, $\varnothing$ 50mm
Mounting quantity per time	2
Moulding temperature	50÷300
Mounting time	<15min
Heating time to 100°C	8min
Heating time to 140°C	10min
Heating time to 180°C	12min
Heating power	1000W
Cooling time	3÷5min
Pressure	0÷72MPa
Noise level	65dB
Operating temperature	0÷30°C
Humidity non-condensing	<85%
Water inlet	8mm
Water outlet	13mm
Power supply	230V $\pm$ 10% - 50/60Hz - 1Ph
Dimensions	630x520x520mm
Weight	120

## Standard accessories

Two molds selected from:  $\varnothing$  25mm,  $\varnothing$  30mm,  $\varnothing$  40mm,  $\varnothing$  50mm; mounting powder; filler; inlet water pipe; outlet water pipe; hose clamp