


DIAMONDCUT-25

Features

- > Bench-top automatic diamond cutting machine, suitable for hard materials, especially suggested for artificial crystals
- > Equipped with four different vices, allows the best cutting angle
- > Easy to use with a compact body and safety protection
- > High level of positioning accuracy
- > Built-in cooling system with tank



Model	DIAMONDCUT-25	
Cutting capacity	Max. \varnothing 25mm	
Cutting disc diameter	100x0.5x12.7mm	
X axis travel	12mm, step 0.1mm	
Noise level	<65dB	
Working chamber LxW	196x104mm	
Speed	50÷800rpm	
Operating temperature	5÷30°C	
Humidity non-condensing	≤85%	
Motor power	100W	
Power supply	230V ±10% - 50/60Hz - 1Ph	
Tank capacity	0.5l	
Machine	Dimensions	350x380x290mm
	Weight	25kg

Standard accessories

Diamond cutting disc \varnothing 100mm; Spanners